


# Cabinet Dry Glove Box

- ✓ **Close system for control Environment.**
- ✓ **Stainless Steel SUS 304 Workstations.**
- ✓ **Chemical Resistant Easy to Cleaning.**
- ✓ **User protection.**

## Application

- Cabinet Close system testing for Laboratories of pharmaceuticals, Biological, Chemical.
- Weighing samples of Drug powders.
- Operation of various instruments requiring environmental control, such as Automatic titration, Volumetric/Culometric Karl Fischer Titration.



ធុនស្រទាប់ HEPA/Carbon FILTER






UVC Lighting 15 W



DRY GLOVES PIERCAN

Options.





## รายละเอียดคุณลักษณะ

1. ตู้ปฏิบัติงานพร้อมถุงมือ ควบคุมสภาวะแวดล้อมการทำงาน เหมาะสำหรับเครื่องมือ อุปกรณ์ สารเคมี และเตรียมตัวอย่างอันตราย
2. ทำจากอะครีลิคใส ความหนา 10 มิลลิเมตร
3. พื้นที่ปฏิบัติงาน ทำจากสแตนเลส 304 ขนาด 60 x 55 ซม. ทนสารเคมี และง่ายต่อการดูแลรักษา
4. หลอดอุลตราไวโอเลต (UVC) ขนาด 15 วัตต์ จำนวน 1 หลอด (Options.)
5. หลอดแสงสว่างฟลูออเรสเซนต์ ขนาด 15 วัตต์ จำนวน 1 หลอด (Options.)
6. ช่องสวมถุงมือยาว ด้านหน้า พร้อมตัวล็อกถุงมือ
7. ปลั๊กไฟ 2 ช่อง สำหรับต่อฟองอุปกรณ์เครื่องมือ
8. ใช้ไฟฟ้า 220 โวลท์ 50 เฮิรตซ์
9. รับประกันคุณภาพ 1 ปี

MODEL	Width (W)	Depth (D)	Height (H)	Area work station.
<b>Cabinet-G60B-30</b>	60 cm.	55 cm.	65 cm.	60 x 55 cm.
<b>Cabinet-G100B-30</b>	100 cm.	55 cm.	65 cm.	100 x 55 cm.